

CDK5RAP3 Rabbit mAb

Catalog No: #57209

Package Size: #57209-1 50ul #57209-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	CDK5RAP3 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IHC IP
Species Reactivity	Human,Mouse
Immunogen Description	A synthetic peptide of human CDK5RAP3
Other Names	C53; IC53; LZAP; HSF-27; MST016; PP1553; OK/SW-cl.114
Accession No.	Swiss-Prot:Q96JB5 Gene ID:80279
Uniprot	Q96JB5
GeneID	80279;
Calculated MW	57 kDa
SDS-PAGE MW	70 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IP:1/20WB: 1/500-1/1000IHC: 1/50-1/100

Images



Western blot analysis of CDK5RAP3 in K562, Hela lysates using CDK5RAP3 antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using CDK5RAP3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

Note: This product is for in vitro research use only