

ADH5 Rabbit mAb

Catalog No: #57214

Package Size: #57214-1 50ul #57214-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	ADH5 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IHC IP
Species Reactivity	Human
Immunogen Description	A synthetic peptide of human ADH5
Other Names	FDH; ADHX; ADH-3; FALDH; GSNOR; GSH-FDH; HEL-S-60p
Accession No.	Swiss-Prot:P11766 Gene ID:128
Uniprot	P11766
GeneID	128;
Calculated MW	40 kDa
SDS-PAGE MW	40 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

IP: 1/20WB: 1/500-1/1000IHC: 1/50-1/100

Images



Western blot analysis of ADH5 in MCF-7, HepG2, K562 lysates using ADH5 antibody.

Note: This product is for in vitro research use only