

ALKBH4 Rabbit mAb

Catalog No: #57217

Package Size: #57217-1 50ul #57217-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ALKBH4 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IP
Species Reactivity	Human
Immunogen Description	Recombinant protein of human ALKBH4
Other Names	ABH4
Accession No.	Swiss-Prot:Q9NXW9 Gene ID:54784
Uniprot	Q9NXW9
GeneID	54784;
Calculated MW	34 kDa
SDS-PAGE MW	34 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

IP:1/20WB: 1/500-1/1000

Images



Western blot analysis of ALKBH4 in K562, HeLa lysates using ALKBH4 antibody.

Note: This product is for in vitro research use only