

## AGPS Rabbit mAb

Catalog No: #57218

Package Size: #57218-1 50ul #57218-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	AGPS Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IHC IP
Species Reactivity	Human
Immunogen Description	Recombinant protein of human Alkyl-DHAP synthase
Other Names	ADAS; ADPS; RCDP3; ADAP-S; ADHAPS; ALDHPSY
Accession No.	Swiss-Prot:O00116 Gene ID:8540
Uniprot	O00116
GeneID	8540;
Calculated MW	73 kDa
SDS-PAGE MW	73 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


## Application Details

IP:1/20WB: 1/500-1/1000IHC: 1/50-1/100

## Images



Western blot analysis of AlkylDHAP synthase in K562 lysates using AGPS antibody.



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using AlkyIDHAP synthase antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

---

Note: This product is for in vitro research use only