

Angiotensinogen Rabbit mAb

Catalog No: #57219

Package Size: #57219-1 50ul #57219-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

| | |
|-----------------------|---|
| Product Name | Angiotensinogen Rabbit mAb |
| Host Species | Rabbit |
| Clonality | Monoclonal |
| Isotype | Rabbit IgG |
| Applications | WB IHC IP |
| Species Reactivity | Mouse |
| Immunogen Description | Recombinant protein of mouse Angiotensinogen |
| Other Names | Angiotensinogen; Serpin A8; AGT; SERPINA8; Des-Asp[1]-angiotensin II; Serine (or cysteine) proteinase inhibitor |
| Accession No. | Swiss-Prot:P11859 Gene ID:11606 |
| Uniprot | P11859 |
| GeneID | 11606; |
| Calculated MW | 52 kDa |
| SDS-PAGE MW | 52 kDa |
| Formulation | 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA |
| Storage | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. |

Application Details

IP:1/20WB: 1/500-1/1000IHC: 1/50-1/100

Images



Western blot analysis of Angiotensinogen in mouse liver lysates using Angiotensinogen antibody.

Note: This product is for in vitro research use only