

## APPBP1 Rabbit mAb

Catalog No: #57221

Package Size: #57221-1 50ul #57221-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	APPBP1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IP
Species Reactivity	Human
Immunogen Description	A synthetic peptide of human APPBP1
Other Names	NEDD8 activating enzyme E1 subunit 1; HPP1; ula-1; APPBP1; A-116A10.1
Accession No.	Swiss-Prot:Q13564 Gene ID:8883
Uniprot	Q13564
GeneID	8883;
Calculated MW	60 kDa
SDS-PAGE MW	60 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

## Application Details

IP:1/20WB: 1/500-1/1000

## Images



Western blot analysis of APPBP1 in SKOV-3 lysates using APPBP1 antibody.

Note: This product is for in vitro research use only