

Ataxin 1 Rabbit mAb

Catalog No: #57224

Package Size: #57224-1 50ul #57224-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Ataxin 1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IP
Species Reactivity	Mouse
Immunogen Description	Recombinant protein of human Ataxin 1
Other Names	ATX1; SCA1; D6S504E
Accession No.	Swiss-Prot:P54253 Gene ID:6310
Uniprot	P54253
GeneID	6310;
Calculated MW	87 kDa
SDS-PAGE MW	105 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

IP:1/20WB: 1/500-1/1000

Images



Western blot analysis of Ataxin 1 in 3T3 lysates using Ataxin 1 antibody.

Note: This product is for in vitro research use only