

B7H4 Rabbit mAb

Catalog No: #57227

Package Size: #57227-1 50ul #57227-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	B7H4 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IHC ICC
Species Reactivity	Human
Immunogen Description	Recombinant protein of human B7H4
Other Names	B7X; B7H4; B7S1; B7-H4; B7h.5; VCTN1; PRO1291
Accession No.	Swiss-Prot:Q7Z7D3 Gene ID:79679
Uniprot	Q7Z7D3
GeneID	79679;
Calculated MW	31 kDa
SDS-PAGE MW	70 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB: 1/500-1/1000 IHC: 1/50-1/100 ICC: 1/50-1/200

Images



Western blot analysis of B7H4 in MDA-MB-468, T47D lysates using B7H4 antibody.

Note: This product is for in vitro research use only