

## BRG1 Rabbit mAb

Catalog No: #57230

Package Size: #57230-1 50ul #57230-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	BRG1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IHC IP ICC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	A synthetic peptide of mouse BRG1
Other Names	Brg1; BAF190A; SW1/SNF; HP1-BP72; SNF2beta; b2b692Clo; b2b508.1Clo
Accession No.	Swiss-Prot:Q3TKT4 Gene ID:20586
Uniprot	Q3TKT4
GeneID	20586;
Calculated MW	181 kDa
SDS-PAGE MW	220 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


## Application Details

IP: 1/20WB: 1/500-1/1000IHC: 1/50-1/100ICC: 1/50-1/200

## Images



Western blot analysis of BRG1 in 3T3, Raw264.7, C6 lysates using BRG1 antibody.



Western blot analysis of BRG1 in 293T, HeLa, HepG2, A549 lysates using BRG1 antibody.

Note: This product is for in vitro research use only