

Calbindin Rabbit mAb

Catalog No: #57232

Package Size: #57232-1 50ul #57232-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Calbindin Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IHC IP
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant protein of mouse Calbindin
Other Names	Calbindin
Accession No.	Swiss-Prot:P12658 Gene ID:12307
Uniprot	P12658
GeneID	12307;
Calculated MW	30 kDa
SDS-PAGE MW	30 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IP: 1/20WB: 1/500-1/1000IHC: 1/50-1/100

Images



Western blot analysis of Calbindin in rat Brain lysates using Calbindin antibody.



Immunohistochemistry analysis of paraffin-embedded mouse kidney using Calbindin antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

Note: This product is for in vitro research use only