

Caspase 11 Rabbit mAb

Catalog No: #57234

Package Size: #57234-1 50ul #57234-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Caspase 11 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IP
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant protein of mouse Caspase-11
Other Names	Casp1; ich-3; CASP-4; Casp11; CASP-11
Accession No.	Swiss-Prot:P70343 Gene ID:12363
Uniprot	P70343
GeneID	12363;
Calculated MW	43 kDa
SDS-PAGE MW	43,38,25 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

IP:1/20WB: 1/500-1/1000

Images



Western blot analysis of Caspase 11 in C6, 3T3, Hela lysates using Caspase 11 antibody.

Note: This product is for in vitro research use only