

Caspase 9 Rabbit mAb

Catalog No: #57235

Package Size: #57235-1 50ul #57235-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Caspase 9 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IP
Species Reactivity	Mouse,Rat
Immunogen Description	Recombinant protein of mouse Caspase-9
Other Names	CASP 9; Caspase-9; CASP-9; Apoptotic protease Mch-6; Apoptotic protease-activating factor 3; APAF-3; ICE-like apoptotic protease 6; ICE-LAP6; Caspase-9 subunit p35; Caspase-9 subunit p10; CASP9; MCH6
Accession No.	Swiss-Prot:Q8C3Q9 Gene ID:12371
Uniprot	Q8C3Q9
GeneID	12371;
Calculated MW	50 kDa
SDS-PAGE MW	50 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IP:1/20WB: 1/500-1/1000

Images



Western blot analysis of Caspase9 in 3T3 lysates using Caspase 9 antibody.



Western blot analysis of Caspase9 in C6 lysates using Caspase9 antibody.

Note: This product is for in vitro research use only