

CD276 Rabbit mAb

Catalog No: #57240

Package Size: #57240-1 50ul #57240-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CD276 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IHC IP
Species Reactivity	Human
Immunogen Description	Recombinant protein of mouse CD276
Other Names	B7h3; B7RP-2; AU016588; 6030411F23Rik
Accession No.	Swiss-Prot:Q8VE98 Gene ID:102657
Uniprot	Q8VE98
GeneID	102657;
Calculated MW	34 kDa
SDS-PAGE MW	90-110 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

IP: 1/20WB: 1/500-1/1000IHC: 1/50-1/100

Images



Western blot analysis of CD276 in HepG2 lysates using CD276 antibody.

Note: This product is for in vitro research use only