

## CDKN2A/p16INK4a Rabbit mAb

Catalog No: #57241

Package Size: #57241-1 50ul #57241-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	CDKN2A/p16INK4a Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IP ICC
Species Reactivity	Mouse,Rat
Immunogen Description	Recombinant protein of mouse CDKN2A/p16INK4a
Other Names	Cdkn2a;P16ink4a;Cyclin-dependent kinase 4 inhibitor A (CDK4I);p16-INK4a ;p16-INK4
Accession No.	Swiss-Prot:P51480 Gene ID:12578
Uniprot	P51480
GeneID	12578;
Calculated MW	18 kDa
SDS-PAGE MW	18 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

## Application Details


IP: 1/20WB: 1/500-1/1000ICC: 1/50-1/200

## Images



Western blot analysis of CDKN2A/p16INK4a in 3T3 lysates using CDKN2A antibody.

Western blot analysis of CDKN2A/p16INK4a in C6 lysates using CDKN2A/p16INK4a antibody.



Note: This product is for in vitro research use only