

CEL2 Rabbit mAb

Catalog No: #57242

Package Size: #57242-1 50ul #57242-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	CEL2 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IHC IP ICC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant protein of human CELF-2
Other Names	ETR3; ETR-3; NAPOR; CELF-2; CUGBP2; BRUNOL3; CUG-BP2
Accession No.	Swiss-Prot:O95319 Gene ID:10659
Uniprot	O95319
GeneID	10659;
Calculated MW	54 kDa
SDS-PAGE MW	54 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IP: 1/20WB: 1/500-1/1000IHC: 1/50-1/100ICC: 1/50-1/200

Images



Western blot analysis of CELF2 in mouse brain, rat hippocampus lysates using CELF2 antibody.



Western blot analysis of CELF2 in HepG2, 293, HeLa, K562 lysates using CELF2 antibody.

Note: This product is for in vitro research use only