

D4 GDI Rabbit mAb

Catalog No: #57251

Package Size: #57251-1 50ul #57251-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	D4 GDI Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IHC ICC
Species Reactivity	Human
Immunogen Description	A synthetic peptide of human D4 GDI
Other Names	Rho GDP dissociation inhibitor (GDI) beta; D4; GDIA2; GDID4; LYGDI; Ly-GDI; RAP1GN1; RhoGDI2
Accession No.	Swiss-Prot:P52566 Gene ID:397
Uniprot	P52566
GeneID	397;
Calculated MW	23 kDa
SDS-PAGE MW	23 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB: 1/500-1/1000 IHC: 1/50-1/100 ICC: 1/50-1/200

Images



Western blot analysis of D4 GDI in Jurkat, U937 lysates using D4 GDI antibody.

Note: This product is for in vitro research use only