

DAP1 Rabbit mAb

Catalog No: #57252

Package Size: #57252-1 50ul #57252-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	DAP1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IHC IP ICC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	A synthetic peptide of human DAP1
Other Names	Death-associated protein 1; DAP-1
Accession No.	Swiss-Prot:P51397 Gene ID:1611
Uniprot	P51397
GeneID	1611;
Calculated MW	11 kDa
SDS-PAGE MW	11 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

IP:1/20WB: 1/500-1/1000IHC: 1/50-1/100ICC: 1/50-1/200

Images



Western blot analysis of DAP1 in MCF-7, C6, 3T3 lysates using DAP1 antibody.

Note: This product is for in vitro research use only