

## DBC 1 Rabbit mAb

Catalog No: #57253

Package Size: #57253-1 50ul #57253-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	DBC 1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IHC IP ICC
Species Reactivity	Human
Immunogen Description	A synthetic peptide of human DBC-1
Other Names	DBC1; DBC-1; NET35; p30DBC; p30 DBC; KIAA1967
Accession No.	Swiss-Prot:Q8N163 Gene ID:57805
Uniprot	Q8N163
GeneID	57805;
Calculated MW	103 kDa
SDS-PAGE MW	130 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

## Application Details

IP: 1/20WB: 1/500-1/1000IHC: 1/50-1/100ICC: 1/50-1/200

## Images



Western blot analysis of DBC1 in 293, Jurkat lysates using DBC 1 antibody.

Note: This product is for in vitro research use only