

DNAJC8 Rabbit mAb

Catalog No: #57255

Package Size: #57255-1 50ul #57255-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	DNAJC8 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IHC IP ICC
Species Reactivity	Human,Mouse,Rat,Monkey
Immunogen Description	A synthetic peptide of human DjC8
Other Names	SPF31; HSPC331
Accession No.	Swiss-Prot:O75937 Gene ID:22826
Uniprot	O75937
GeneID	22826;
Calculated MW	30 kDa
SDS-PAGE MW	30 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IP: 1/20WB: 1/500-1/1000IHC: 1/50-1/100ICC: 1/50-1/200

Images



Western blot analysis of DjC8 in mouse brain, rat lung lysates using DNAJC8 antibody.



Western blot analysis of DjC8 in Jurkat, HeLa, COS7, Raji lysates using DjC8 antibody.

Note: This product is for in vitro research use only