

## DLX4 Rabbit mAb

Catalog No: #57256



Package Size: #57256-1 50ul #57256-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	DLX4 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IP ICC
Species Reactivity	Human,Rat
Immunogen Description	A synthetic peptide of human Dlx4
Other Names	DLX4; BP1; DLX7; DLX8; DLX9; Homeobox protein DLX-4; Beta protein 1; Homeobox protein DLX-7; Homeobox protein DLX-8
Accession No.	Swiss-Prot:Q92988 Gene ID:1748
Uniprot	Q92988
GeneID	1748;
Calculated MW	26 kDa
SDS-PAGE MW	35 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

## Application Details

IP:1/20WB: 1/500-1/1000ICC: 1/50-1/200

Note: This product is for in vitro research use only