

## FAM98B Rabbit mAb

Catalog No: #57260

Package Size: #57260-1 50ul #57260-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	FAM98B Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IHC ICC
Species Reactivity	Human
Immunogen Description	Recombinant protein of human FAM98B
Other Names	Protein FAM98B
Accession No.	Swiss-Prot:Q52LJ0 Gene ID:283742
Uniprot	Q52LJ0
GeneID	283742;
Calculated MW	46 kDa
SDS-PAGE MW	37 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


## Application Details

WB: 1/500-1/1000 IHC: 1/50-1/100 ICC: 1/50-1/200


## Images



Western blot analysis of FAM98B in K562 lysates using FAM98B antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using FAM98B antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of FAM98B (green) in HeLa using FAM98B antibody, and DAPI (blue)

---

Note: This product is for in vitro research use only