

FCGR1A Rabbit mAb

Catalog No: #57262

Package Size: #57262-1 50ul #57262-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	FCGR1A Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IP ICC
Species Reactivity	Human
Immunogen Description	A synthetic peptide of human FCGR1A
Other Names	Fc-gamma RI; FcRI; Fc-gamma RIA; FcgammaRIa; CD64
Accession No.	Swiss-Prot:P12314 Gene ID:2209
Uniprot	P12314
GeneID	2209;
Calculated MW	43 kDa
SDS-PAGE MW	43-80 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

IP: 1/20WB: 1/500-1/1000ICC: 1/50-1/200

Images



Western blot analysis of FCGR1A in THP-1 lysates using FCGR1A antibody.

Note: This product is for in vitro research use only