

Fibulin 1 Rabbit mAb

Catalog No: #57263

Package Size: #57263-1 50ul #57263-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Fibulin 1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IP
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant protein of human Fibulin 1
Other Names	FBLN; FIBL1
Accession No.	Swiss-Prot:P23142 Gene ID:2192
Uniprot	P23142
GeneID	2192;
Calculated MW	77 kDa
SDS-PAGE MW	77 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

IP:1/20WB: 1/500-1/1000

Images



Western blot analysis of Fibulin 1 in Human placenta, A549, 3T3 lysates using Fibulin 1 antibody.

Note: This product is for in vitro research use only