

FNTA Rabbit mAb

Catalog No: #57264

Package Size: #57264-1 50ul #57264-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	FNTA Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IHC ICC
Species Reactivity	Human,Rat
Immunogen Description	A synthetic peptide of human FNTA
Other Names	FPTA; PTAR2; PGGT1A
Accession No.	Swiss-Prot:P49354 Gene ID:2339
Uniprot	P49354
GeneID	2339;
Calculated MW	44 kDa
SDS-PAGE MW	44 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB: 1/500-1/1000 IHC: 1/50-1/100 ICC: 1/50-1/200

Images



Western blot analysis of FNTA in Jurkat, PC-12, HepG2 lysates using FNTA antibody.

Note: This product is for in vitro research use only