

FOXC1 Rabbit mAb

Catalog No: #57265

Package Size: #57265-1 50ul #57265-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	FOXC1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IP
Species Reactivity	Mouse
Immunogen Description	Recombinant protein of human FOXC1
Other Names	ARA; IGDA; IHG1; FKHL7; IRID1; RIEG3; FREAC3; FREAC-3
Accession No.	Swiss-Prot:Q12948 Gene ID:2296
Uniprot	Q12948
GeneID	2296;
Calculated MW	57 kDa
SDS-PAGE MW	75 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

IP:1/20WB: 1/500-1/1000

Images



Western blot analysis of FOXC1 in 3T3 lysates using FOXC1 antibody.

Note: This product is for in vitro research use only