

## KAT2A Rabbit mAb

Catalog No: #57266

Package Size: #57266-1 50ul #57266-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	KAT2A Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IP
Species Reactivity	Human
Immunogen Description	Recombinant protein of human GCN5
Other Names	GCN5; hGCN5; GCN5L2; PCAF-b
Accession No.	Swiss-Prot:Q92830 Gene ID:2648
Uniprot	Q92830
GeneID	2648;
Calculated MW	94 kDa
SDS-PAGE MW	94 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

## Application Details

IP:1/20WB: 1/500-1/1000

## Images



Western blot analysis of KAT2A / GCN5 in K562 lysates using KAT2A antibody.

Note: This product is for in vitro research use only