

GFPT2 Rabbit mAb

Catalog No: #57267

Package Size: #57267-1 50ul #57267-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	GFPT2 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IP ICC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant protein of human GFPT2
Other Names	GFAT; GFAT2; GFAT 2
Accession No.	Swiss-Prot:O94808 Gene ID:9945
Uniprot	O94808
GeneID	9945;
Calculated MW	77 kDa
SDS-PAGE MW	77 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

IP: 1/20WB: 1/500-1/1000ICC: 1/50-1/200


Images



Western blot analysis of GFPT2 in rat spleen lysates using GFPT2 antibody.



Western blot analysis of GFPT2 in 3T3, HeLa lysates using GFPT2 antibody



Immunocytochemistry analysis of GFPT2 (green) in HT-1080 using GFPT2 antibody, and DAPI(blue).

Note: This product is for in vitro research use only