

## HNRNPA1L2 Rabbit mAb

Catalog No: #57269

Package Size: #57269-1 50ul #57269-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	HNRNPA1L2 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IHC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	A synthetic peptide of human Heterogeneous Nuclear Ribonucleoprotein (A1-like)
Other Names	hnRNP core protein A1-like 2
Accession No.	Swiss-Prot:Q32P51 Gene ID:144983
Uniprot	Q32P51
GeneID	144983;
Calculated MW	34 kDa
SDS-PAGE MW	34 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

## Application Details

WB: 1/500-1/1000IHC: 1/50-1/100

## Images



Western blot analysis of Heterogeneous Nuclear Ribonucleoprotein (A1like) in K562, C6, 3T3, HeLa lysates using HNRNPA1L2 antibody.

Note: This product is for in vitro research use only