

HLA C Rabbit mAb

Catalog No: #57272

Package Size: #57272-1 50ul #57272-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	HLA C Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IP
Species Reactivity	Human,Mouse
Immunogen Description	A synthetic peptide of human HLA-C
Other Names	MHC; HLAC; HLC-C; D6S204; PSORS1; HLA-JY3
Accession No.	Swiss-Prot:P10321 Gene ID:3107
Uniprot	P10321
GeneID	3107;
Calculated MW	41 kDa
SDS-PAGE MW	41 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

IP:1/20WB: 1/500-1/1000

Images



Western blot analysis of HLAC in HeLa, 293T, Jurkat, 3T3, HepG2 lysates using HLA C antibody.

Note: This product is for in vitro research use only