

HuC/HuD Rabbit mAb

Catalog No: #57277

Package Size: #57277-1 50ul #57277-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	HuC/HuD Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IP
Species Reactivity	Rat
Immunogen Description	Recombinant protein of mouse HuC
Other Names	Huc; mHuC; PLE21; 2600009P04Rik
Accession No.	Swiss-Prot:Q60900/P26378 Gene ID:15571/1996
Uniprot	Q60900/P26378
GeneID	15571/1996;
Calculated MW	40 kDa
SDS-PAGE MW	40 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

IP:1/20WB: 1/500-1/1000

Images



Western blot analysis of HuD/HuC in rat Brain lysates using HuC/HuD antibody.

Note: This product is for in vitro research use only