

ICOS Ligand Rabbit mAb

Catalog No: #57278

Package Size: #57278-1 50ul #57278-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	ICOS Ligand Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IP
Species Reactivity	Human
Immunogen Description	Recombinant protein of human ICOS Ligand
Other Names	inducible T-cell co-stimulator ligand; B7H2; GL50; B7-H2; B7RP1; CD275; ICOSL; LICOS; B7RP-1; ICOS-L
Accession No.	Swiss-Prot:O75144 Gene ID:23308
Uniprot	O75144
GeneID	23308;
Calculated MW	33 kDa
SDS-PAGE MW	75 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

IP:1/20WB: 1/500-1/1000

Images



Western blot analysis of ICOS Ligand in K562 lysates using ICOS Ligand antibody.

Note: This product is for in vitro research use only