

kappa Opioid Receptor Rabbit mAb

Catalog No: #57282

Package Size: #57282-1 50ul #57282-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	kappa Opioid Receptor Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IP
Species Reactivity	Rat
Immunogen Description	A synthetic peptide of mouse Kappa Opioid Receptor
Other Names	KOR; R21; KOR-1; MSL-1; Oprk2; K-OR-1
Accession No.	Swiss-Prot:P33534 Gene ID:18387
Uniprot	P33534
GeneID	18387;
Calculated MW	43 kDa
SDS-PAGE MW	60 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

IP:1/20WB: 1/500-1/1000

Images



Western blot analysis of Kappa Opioid Receptor in rat Brain lysates using kappa Opioid Receptor antibody.

Note: This product is for in vitro research use only