

## LPCAT1 Rabbit mAb

Catalog No: #57284

Package Size: #57284-1 50ul #57284-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	LPCAT1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IHC IP ICC
Species Reactivity	Human
Immunogen Description	Recombinant protein of human LPCAT1
Other Names	AYTL2; lpcat; AGPAT9; PFAAP3; AGPAT10; LPCAT-1; lysoPAFAT
Accession No.	Swiss-Prot:Q8NF37 Gene ID:79888
Uniprot	Q8NF37
GeneID	79888;
Calculated MW	59 kDa
SDS-PAGE MW	29 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

## Application Details

IP: 1/20WB: 1/500-1/1000IHC: 1/50-1/100ICC: 1/50-1/200

## Images



Western blot analysis of LPCAT1 in A549 lysates using LPCAT1 antibody.

Note: This product is for in vitro research use only