

MAPK6 Rabbit mAb

Catalog No: #57286



Package Size: #57286-1 50ul #57286-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	MAPK6 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IHC IP
Species Reactivity	Human,Mouse
Immunogen Description	A synthetic peptide of human MAPK6
Other Names	ERK3; PRKM6; p97MAPK; HsT17250; ERK-3
Accession No.	Swiss-Prot:Q16659 Gene ID:5597
Uniprot	Q16659
GeneID	5597;
Calculated MW	83 kDa
SDS-PAGE MW	105 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

IP:1/20WB: 1/500-1/1000IHC: 1/50-1/100

Note: This product is for in vitro research use only