

MMP9 Rabbit mAb

Catalog No: #57292

Package Size: #57292-1 50ul #57292-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	MMP9 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IHC IP
Species Reactivity	Rat
Immunogen Description	Recombinant protein of mouse MMP9
Other Names	92 kDa gelatinase; 92 kDa type IV collagenase; Gelatinase B; GELB
Accession No.	Swiss-Prot:P41245 Gene ID:17395
Uniprot	P41245
GeneID	17395;
Calculated MW	81 kDa
SDS-PAGE MW	92 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

IP: 1/20WB: 1/500-1/1000IHC: 1/50-1/100

Images



Western blot analysis of MMP9 in rat lung lysates using MMP9 antibody.

Note: This product is for in vitro research use only