

NCE2 Rabbit mAb

Catalog No: #57297

Package Size: #57297-1 50ul #57297-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	NCE2 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IHC ICC
Species Reactivity	Human
Immunogen Description	Recombinant protein of human NCE2
Other Names	NCE2; NEDD8 carrier protein UBE2F; NEDD8 protein ligase UBE2F; NEDD8-conjugating enzyme 2
Accession No.	Swiss-Prot:Q969M7 Gene ID:140739
Uniprot	Q969M7
GeneID	140739;
Calculated MW	21 kDa
SDS-PAGE MW	21 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB: 1/500-1/1000 IHC: 1/50-1/100 ICC: 1/50-1/200

Images



Western blot analysis of NCE2 in HeLa, MCF-7 lysates using NCE2 antibody.

Note: This product is for in vitro research use only