

NDRG2 Rabbit mAb

Catalog No: #57298

Package Size: #57298-1 50ul #57298-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	NDRG2 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IHC IP ICC
Species Reactivity	Mouse,Rat
Immunogen Description	A synthetic peptide of human NDRG2
Other Names	SYLD
Accession No.	Swiss-Prot:Q9UN36 Gene ID:57447
Uniprot	Q9UN36
GeneID	57447;
Calculated MW	41 kDa
SDS-PAGE MW	41 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details


IP:1/20WB: 1/500-1/1000IHC: 1/50-1/100ICC: 1/50-1/200

Images



Western blot analysis of NDRG2 in rat brain lysates using NDRG2 antibody.

Western blot analysis of NDRG2 in mouse brain lysates using NDRG2 antibody.



Note: This product is for in vitro research use only