

NDUFAF4 Rabbit mAb

Catalog No: #57299



Package Size: #57299-1 50ul #57299-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	NDUFAF4 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IP ICC
Species Reactivity	Human
Immunogen Description	Recombinant protein of human NDUFAF4
Other Names	My013; C6orf66; HRPAP20; HSPC125; bA22L21.1
Accession No.	Swiss-Prot:Q9P032 Gene ID:29078
Uniprot	Q9P032
GeneID	29078;
Calculated MW	20 kDa
SDS-PAGE MW	20 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

IP:1/20WB: 1/500-1/1000ICC: 1/50-1/200

Note: This product is for in vitro research use only