

NDUFAF5 Rabbit mAb

Catalog No: #57300



Package Size: #57300-1 50ul #57300-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	NDUFAF5 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IHC IP
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant protein of human NDUFAF5
Other Names	C20orf7; dJ842G6.1; bA526K24.2
Accession No.	Swiss-Prot:Q5TEU4 Gene ID:79133
Uniprot	Q5TEU4
GeneID	79133;
Calculated MW	39 kDa
SDS-PAGE MW	39 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

IP:1/20WB: 1/500-1/1000IHC: 1/50-1/100

Note: This product is for in vitro research use only