

## Neurogranin Rabbit mAb

Catalog No: #57304

Package Size: #57304-1 50ul #57304-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

|                       |  |
|-----------------------|--|
| Product Name          | Neurogranin Rabbit mAb   |
| Host Species          | Rabbit   |
| Clonality             | Monoclonal   |
| Isotype               | Rabbit IgG   |
| Applications          | WB IHC IP ICC  |
| Species Reactivity    | Human,Mouse,Rat  |
| Immunogen Description | Recombinant protein of human Neurogranin   |
| Other Names           | neurogranin (protein kinase C substrate; RC3); RC3; hng                                  |
| Accession No.         | Swiss-Prot:Q92686 Gene ID:4900   |
| Uniprot               | Q92686   |
| GeneID                | 4900;  |
| Calculated MW         | 8 kDa  |
| SDS-PAGE MW           | 14-16 kDa  |
| Formulation           | 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA    |
| Storage               | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. |

## Application Details

IP: 1/20WB: 1/500-1/1000IHC: 1/50-1/100ICC: 1/50-1/200

## Images



Western blot analysis of Neurogranin in U87-MG, PC-3 lysates using Neurogranin antibody.

Note: This product is for in vitro research use only