

NIT1 Rabbit mAb

Catalog No: #57305

Package Size: #57305-1 50ul #57305-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	NIT1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IP
Species Reactivity	Human,Mouse,Rat
Immunogen Description	A synthetic peptide of human NIT1
Other Names	Nitrilase homolog 1
Accession No.	Swiss-Prot:Q86X76 Gene ID:4817
Uniprot	Q86X76
GeneID	4817;
Calculated MW	36 kDa
SDS-PAGE MW	36 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

IP:1/20WB: 1/500-1/1000

Images



Western blot analysis of NIT1 in HeLa, 293T, Jurkat, 3T3 lysates using NIT1 antibody.

Note: This product is for in vitro research use only