

Olfactory Marker Protein Rabbit mAb

Catalog No: #57306

Package Size: #57306-1 50ul #57306-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Olfactory Marker Protein Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IP
Species Reactivity	Mouse,Rat
Immunogen Description	Recombinant protein of mouse Olfactory Marker Protein
Other Names	Olfactory marker protein
Accession No.	Swiss-Prot:Q64288 Gene ID:18378
Uniprot	Q64288
GeneID	18378;
Calculated MW	19 kDa
SDS-PAGE MW	19 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

IP:1/20WB: 1/500-1/1000

Images



Western blot analysis of Olfactory Marker Protein in rat Brain, MCF-7 lysates using Olfactory Marker Protein antibody.

Note: This product is for in vitro research use only