

## p53 Rabbit mAb

Catalog No: #57308

Package Size: #57308-1 50ul #57308-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	p53 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IP ICC
Species Reactivity	Mouse
Immunogen Description	Recombinant protein of mouse p53
Other Names	bb1; bfy; bhy; p44; p53; Tp53
Accession No.	Swiss-Prot:P02340 Gene ID:22059
Uniprot	P02340
GeneID	22059;
Calculated MW	43 kDa
SDS-PAGE MW	53 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

## Application Details

IP: 1/20WB: 1/500-1/1000ICC: 1/50-1/200

## Images



Western blot analysis of p53 in 3T3 lysates using p53 antibody.

Note: This product is for in vitro research use only