

PAIP1 Rabbit mAb

Catalog No: #57309

Package Size: #57309-1 50ul #57309-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	PAIP1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IHC IP ICC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	A synthetic peptide of human PAIP1
Other Names	poly(A) binding protein interacting protein 1
Accession No.	Swiss-Prot:Q9H074 Gene ID:10605
Uniprot	Q9H074
GeneID	10605;
Calculated MW	54 kDa
SDS-PAGE MW	54 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

IP: 1/20WB: 1/500-1/1000IHC: 1/50-1/100ICC: 1/50-1/200

Images



Western blot analysis of PAIP1 in HeLa, 293T, mouse brain lysates using PAIP1 antibody.

Note: This product is for in vitro research use only