

Pellino 1 Rabbit mAb

Catalog No: #57312

Package Size: #57312-1 50ul #57312-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Pellino 1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IP
Species Reactivity	Human,Rat
Immunogen Description	Recombinant protein of human Pellino 1
Other Names	Pellino-related intracellular-signaling molecule; RING-type E3 ubiquitin transferase pellino homolog 1
Accession No.	Swiss-Prot:Q96FA3 Gene ID:57162
Uniprot	Q96FA3
GeneID	57162;
Calculated MW	46 kDa
SDS-PAGE MW	46 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

IP:1/20WB: 1/500-1/1000

Images



Western blot analysis of Pellino 1 in rat Brain, Jurkat lysates using Pellino 1 antibody.

Note: This product is for in vitro research use only