

## PERK Rabbit mAb

Catalog No: #57313

Package Size: #57313-1 50ul #57313-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	PERK Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant protein of human PERK
Other Names	EIF2AK3; PEK; PERK; Eukaryotic translation initiation factor 2-alpha kinase 3; PRKR-like endoplasmic reticulum kinase; Pancreatic eIF2-alpha kinase; HsPEK
Accession No.	Swiss-Prot:Q9NZJ5 Gene ID:9451
Uniprot	Q9NZJ5
GeneID	9451;
Calculated MW	125 kDa
SDS-PAGE MW	140 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

## Application Details

WB: 1/500-1/1000

## Images



Western blot analysis of PERK in K562, C6, 3T3, HeLa lysates using PERK antibody.

Note: This product is for in vitro research use only