

## PFKP Rabbit mAb

Catalog No: #57314

Package Size: #57314-1 50ul #57314-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	PFKP Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IP ICC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant protein of human PFKP
Other Names	PFKP; PFKF; 6-phosphofructokinase type C; 6-phosphofructokinase; platelet type; Phosphofructo-1-kinase isozyme C; PFK-C; Phosphofructokinase 1; Phosphohexokinase
Accession No.	Swiss-Prot:Q01813 Gene ID:5214
Uniprot	Q01813
GeneID	5214;
Calculated MW	86 kDa
SDS-PAGE MW	86 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

## Application Details


IP:1/20WB: 1/500-1/1000ICC: 1/50-1/200

## Images



Western blot analysis of PFKP in C6, 3T3, HeLa lysates using PFKP antibody.

Immunocytochemistry analysis of PFKP(green) in Hela using PFKP antibody, and DAPI(blue)



---

Note: This product is for in vitro research use only