

PHLDA1 Rabbit mAb

Catalog No: #57315

Package Size: #57315-1 50ul #57315-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	PHLDA1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IHC ICC
Species Reactivity	Human
Immunogen Description	A synthetic peptide of human PHLDA1
Other Names	PQ-rich protein; PQR protein; Proline- and histidine-rich protein; T-cell death-associated gene 51 protein
Accession No.	Swiss-Prot:Q8WV24 Gene ID:22822
Uniprot	Q8WV24
GeneID	22822;
Calculated MW	45 kDa
SDS-PAGE MW	45 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB: 1/500-1/1000 IHC: 1/50-1/100 ICC: 1/50-1/200

Images



Western blot analysis of PHLDA1 in A375 lysates using PHLDA1 antibody.

Note: This product is for in vitro research use only