

PI3 Kinase p110 gamma Rabbit mAb

Catalog No: #57316

Package Size: #57316-1 50ul #57316-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	PI3 Kinase p110 gamma Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IP
Species Reactivity	Human
Immunogen Description	A synthetic peptide of human PI3K-gamma
Other Names	PI3K; PIK3; PI3CG; PI3Kgamma; p110gamma; p120-PI3K
Accession No.	Swiss-Prot:P48736 Gene ID:5294
Uniprot	P48736
GeneID	5294;
Calculated MW	126 kDa
SDS-PAGE MW	110 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

IP:1/20WB: 1/500-1/1000

Images



Western blot analysis of PI3 Kinase p110 gamma in K562 lysates using PI3 Kinase p110 gamma antibody.

Note: This product is for in vitro research use only