

## PIP5K1C Rabbit mAb

Catalog No: #57317

Package Size: #57317-1 50ul #57317-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	PIP5K1C Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Applications	WB IP
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant protein of mouse PIP5KI gamma
Other Names	AI115456; AI835305; Pip5kIgamma
Accession No.	Swiss-Prot:O70161 Gene ID:18717
Uniprot	O70161
GeneID	18717;
Calculated MW	72 kDa
SDS-PAGE MW	90,87 kDa
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


## Application Details

IP:1/20WB: 1/500-1/1000

## Images



Western blot analysis of PIP5K1 gamma in mouse brain lysates using PIP5K1C antibody.



Western blot analysis of PIP5K1 gamma in MCF-7 lysates using PIP5K1 gamma antibody.

---

Note: This product is for in vitro research use only